


LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM



U.S. RETAIL PRODUCT CATALOG AND PRICE LIST*


LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

System Overview and Abutment Selection Criteria

Each LOCATOR R-Tx Attachment System is delivered to you in a convenient, all-in-one double-ended vial that separately holds the Abutment and Processing Components.

Order for your favorite implant systems using this product ordering and reference guide.

LOCATOR R-Tx
Abutment




\$155 PRICE
INCLUDES:

- 1 - Abutment
- 1 - Denture Attachment Housing with Black Processing Insert
- 1 - Standard Zero (**Gray**) Retention Insert
- 1 - Standard Low (**Blue**) Retention Insert
- 1 - Standard Medium (**Pink**) Retention Insert
- 1 - Standard High (**Clear**) Retention Insert
- 1 - Block-Out Spacer

Standard Range Retention Inserts



Zero (**Gray**)
Retention Insert



Low (**Blue**)
Retention Insert




Medium (**Pink**)
Retention Insert



High (**Clear**)
Retention Insert

Limited Range Retention Inserts



Low (**Aqua**)
Retention Insert




Medium (**Magenta**)
Retention Insert



High (**Gold**)
Retention Insert

SELECTING THE PROPER CUFF HEIGHT FOR THE LOCATOR R-Tx ABUTMENT



1. Using a periodontal probe, measure from the top of the implant to the highest point of the gingiva.
2. Select the cuff height of the LOCATOR R-Tx Abutment that corresponds to that measurement. For example, if the measurement is 2mm, select the 2mm cuff. If the measurement is 2.5mm, select the 3mm cuff. An additional 2mm of Abutment height will extend above the gingiva to accommodate for the Denture Attachment Housing and Retention Insert.

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information

LOCATOR R-Tx ABUTMENTS ARE IMPLANT SPECIFIC

TO ORDER, YOU WILL NEED
THE FOLLOWING INFORMATION:

1. Name of the implant system
2. Implant diameter
3. Connection type and platform diameter
4. Tissue measurement (at the highest point)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

ACE SURGICAL

BONE LEVEL CONNECTION (Compatible with infinity Octagon Implants)

2.7mm Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
3.3mm Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06

TISSUE LEVEL CONNECTION (Compatible with infinity Octagon Implants)

4.8mm Regular Platform (RP)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
-----------------------------	--	----------	----------	----------	----------	----------	----------

INTERNAL HEX CONNECTION (Compatible with infinity Internal HEX Implants)

3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06

TRI-CHANNEL CONNECTION (Compatible with infinity TRI-CAM Implants)

3.5mm Diameter Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Diameter Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Diameter Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

EXTERNAL HEX CONNECTION (Compatible with infinity EXTERNAL HEX Implants)

3.5mm Diameter Platform		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Diameter Platform		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm Diameter Platform (Platform Switch Only)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06

AVINENT

INTERNAL HEX CONNECTION (Compatible with Biomimetic Ocean and Coral)

3.5mm I.C (Ocean)	32500-00 (0.5mm)	32500-01	32500-02	32500-03	32500-04	32500-05	32500-06
4.1mm I.C (Ocean & Coral)	32500-00 (0.5mm)	32500-01	32500-02	32500-03	32500-04	32500-05	32500-06

INTERNAL HEX CONNECTION (Compatible with Biomimetic Ocean and Coral)

3.5mm E.C		32505-01	32505-02	32505-03	32505-04	32505-05	32505-06
-----------	--	----------	----------	----------	----------	----------	----------

BIOHORIZONS

INTERNAL HEX CONNECTION (Compatible with Tapered Plus, Tapered Internal and Single-Stage)

3.0mm Diameter Platform		30600-01	30600-02	30600-03	30600-04	30600-05	30600-06
-------------------------	--	----------	----------	----------	----------	----------	----------

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

BIOHORIZONS

3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06

BIOMET 3i NOTE: BIOMET 3i Abutments feature a .048 hex drive

INTERNAL HEX CERTAIN® CONNECTION (Compatible with 3i T3®, NanoTite™ and OSSEOTITE® Implants)

3.4mm Diameter Certain Platform		30300-01	30300-02	30300-03	30300-04	30300-05	30300-06
4.1mm Diameter Certain Platform		30301-01	30301-02	30301-03	30301-04	30301-05	30301-06

EXTERNAL HEX CONNECTION (Compatible with 3i T3®, NanoTite™ and OSSEOTITE® Implants)

3.4mm Diameter External Hex Platform		30305-01	30305-02	30305-03	30305-04	30305-05	30305-06
4.1mm Diameter External Hex Platform		30306-01 (1.2mm)	30306-02	30306-03	30306-04	30306-05	30306-06

BLUE SKY BIO

BONE LEVEL CONNECTION (Compatible with Quattro Implants)

2.8mm Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
3.3mm Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06

TISSUE LEVEL CONNECTION (Compatible with One Stage Implants)

4.8mm Regular Platform (RP)		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06
6.5mm Wide Platform (WP)		30907-01	30907-02	30907-03	30907-04	30907-05	30907-06

DENTSPLY

CONICAL CONNECTION (Compatible with Ankylos® C/X Implants)

2.8mm Ankylos Implant		30715-01	30715-02	30715-03	30715-04	30715-05	30715-06
-----------------------	--	----------	----------	----------	----------	----------	----------

(Compatible with AstraTech Osseospeed EV Implants)

3.6mm Diameter Implant		30700-01	30700-02	30700-03	30700-04	30700-05	30700-06
4.2mm Diameter Implant		30701-01	30701-02	30701-03	30701-04	30701-05	30701-06
4.8mm Diameter Implant		30702-01	30702-02	30702-03	30702-04	30702-05	30702-06

(Compatible with AstraTech Osseospeed TX Implants)

2.9mm Small Platform (Aqua)		30706-01	30706-02	30706-03	30706-04	30706-05	30706-06
3.9mm Large Platform (Lilac)		30707-01	30707-02	30707-03	30707-04	30707-05	30707-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

DENTSPLY

INTERNAL HEX CONNECTION (Compatible with XiVE® Frialit-2 Implants)

3.4mm Diameter Implant		30711-01	30711-02	30711-03	30711-04	30711-05	30711-06
3.8mm Diameter Implant		30712-01	30712-02	30712-03	30712-04	30712-05	30712-06
4.5mm Diameter Implant		30713-01	30713-02	30713-03	30713-04	30713-05	30713-06
5.5mm Diameter Implant		30714-01	30714-02	30714-03	30714-04	30714-05	30714-06

HIOSEN/OSSTEM

EXTERNAL HEX CONNECTION (Compatible with US)

3.5mm Mini Platform (MP) (Yellow)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Regular Platform (RP) (Green)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
5.0mm Wide Platform Switch (WPS) (Blue)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.1mm Wide Platform (WP) (Blue)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06

CONICAL/MORSE TAPE CONNECTION (Compatible with ET and TS)

3.35mm Mini Regular Platform (RP) (Yellow)		31601-01	31601-02	31601-03	31601-04	31601-05	31601-06
--	--	----------	----------	----------	----------	----------	----------

INTERNAL HEX/MORSE TAPER CONNECTION (Compatible with ET and TS)

3.28mm Mini Platform (MP) (Yellow)		31600-01	31600-02	31600-03	31600-04	31600-05	31600-06
------------------------------------	--	----------	----------	----------	----------	----------	----------

INTERNAL OCTAGON/MORSE TAPER CONNECTION (Compatible with SS)

4.8mm Regular Platform (RP) (Green)		31606-01	31606-02	31606-03	31606-04	31606-05	31606-06
-------------------------------------	--	----------	----------	----------	----------	----------	----------

IDS

INTERNAL HEX CONNECTION (Compatible with MegaGen AnyRidge)

3.5mm Small Platform (SP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Regular Platform (RP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Wide Platform (WP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
5.0mm Super Wide Platform		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06

IMPLANT DIRECT

INTERNAL HEX CONNECTION (Compatible with Legacy1, Legacy2, Legacy3 and Legacy4 Implants)

3.0mm Diameter Platform (Blue) (Legacy 2,3,4 Only)		31500-01 (0.75mm)	31500-02	31500-03	31500-04	31500-05	31500-06
3.5mm Diameter Platform (Green)	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform (Purple)	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform (Yellow)		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
4.8mm Diameter Platform		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

IMPLANT DIRECT

TRI-LOBE CONNECTION (Compatible with RePlant and ReActive Implants)

3.5mm Diameter Platform (Purple)	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Diameter Platform (Yellow)	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Diameter Platform (Blue)		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

CONICAL CONNECTION (Compatible with Interactive Implants)

3.0mm Diameter Platform (Purple)		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
3.4mm Diameter Platform (Yellow)		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06

INTERNAL OCTAGON CONNECTION (Compatible with Swish Plus and Swish Tapered Implants)

4.8mm Diameter Platform		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06
-------------------------	--	----------	----------	----------	----------	----------	----------

KEYSTONE DENTAL

TILOBE CONNECTION (Compatible with Genesis and PrimaConnex)

3.0mm Diameter (Green) Platform		31200-01	31200-02	31200-03	31200-04	31200-05	31200-06
---------------------------------	--	----------	----------	----------	----------	----------	----------

(Compatible with PrimaConnex)

3.6mm Diameter (Blue) Platform		31201-01	31201-02	31201-03	31201-04	31201-05	31201-06
--------------------------------	--	----------	----------	----------	----------	----------	----------

(Compatible with PrimaConnex and TILOBEMAXX)

4.3mm Diameter (Purple) Platform		31202-01	31202-02	31202-03	31202-04	31202-05	31202-06
----------------------------------	--	----------	----------	----------	----------	----------	----------

MIS

CONICAL CONNECTION (Compatible with C1 and V3 Implants)

2.75mm Narrow Platform (NP)		32601-01	32601-02	32601-03	32601-04	32601-05	32601-06
3.15mm Standard Platform (SP)		32602-01	32602-02	32602-03	32602-04	32602-05	32602-06

INTERNAL HEX CONNECTION (Compatible with SEVEN and M4 Implants)

3.3mm Narrow Platform (NP)		32605-01	32605-02	32605-03	32605-04	32605-05	32605-06
3.75mm Standard Platform (SP)	32606-00 (0.25mm)	32606-01	32606-02	32606-03	32606-04	32606-05	32606-06
4.5mm Wide Platform (WP)	32607-00 (0.5mm)	32607-01	32607-02	32607-03	32607-04	32607-05	32607-06

NEODENT

CONICAL/CONE MORSE CONNECTION (Compatible with CM Implants)

2.5mm		30715-01	30715-02	30715-03	30715-04	30715-05	30715-06
-------	--	----------	----------	----------	----------	----------	----------

CONICAL/GRAND MORSE CONNECTION (Compatible with GM Implants)

3.0mm		32700-01	32700-02	32700-03	32700-04	32700-05	32700-06
-------	--	----------	----------	----------	----------	----------	----------

NEW

NEW

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

NOBEL BIO CARE

CONICAL CONNECTION (Compatible with NobelActive®, NobelParallel™ CC and NobelReplace® CC Implants)

3.0mm Narrow Platform (NP)		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
3.5mm Regular Platform (RP)		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06
5.1mm Wide Platform (WP)		30508-01	30508-02	30508-03	30508-04	30508-05	30508-06

EXTERNAL HEX CONNECTION (Compatible with Brånemark and NobelSpeedy® Groovy Implants)

3.5mm Narrow Platform (NP)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Regular Platform (RP)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm Wide Platform (WP)		30512-01	30512-02	30512-03	30512-04	30512-05	

TRI-CHANNEL CONNECTION (Compatible with NobelReplace®, ReplaceSelect™ and NobelSpeedy® Replace Implants)

3.5mm Narrow Platform (NP)	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.0mm/4.3mm Regular Platform (RP)	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Wide Platform (WP)		30502-01	30502-02	30502-03	30502-04	30502-05	

OCO BIOMEDICAL

INTERNAL HEX CONNECTION (Compatible with Engage® Implants)

3.0mm Diameter Platform		31500-01 (0.75mm)	31500-02	31500-03	31500-04	31500-05	31500-06
3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06

OSTEOREADY/LITTLE IMPLANT COMPANY

INTERNAL HEX CONNECTION (Compatible with Marc Nevins®, Pamela Ray™ and Performance®)

3.5mm Diameter Platform	32401-00 (0.5mm)	32401-01	32401-02	32401-03	32401-04	32401-05	32401-06
-------------------------	---------------------	----------	----------	----------	----------	----------	----------

SOUTHERN IMPLANTS

TRI-NEX CONNECTION (Compatible with Tri-Nex™ Implants)

3.5mm Tri-Nex Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Tri-Nex Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Tri-Nex Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

EXTERNAL HEX CONNECTION

4.05 Diameter Platform		31703-01	31703-02	31703-03	31703-04	31703-05	31703-06
5.0 Diameter Platform		31704-01	31704-02	31704-03	31704-04	31704-05	31704-06
6.0 Diameter Platform		31705-01	31705-02	31705-03	31705-04	31705-05	31705-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

STRAUMANN

BONE LEVEL (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

3.3mm Narrow Crossfit (NC)		30401-01	30401-02	30401-03	30401-04	30401-05	30401-06
4.1mm/4.8mm Regular Crossfit (RC)		30402-01	30402-02	30402-03	30402-04	30402-05	30402-06

NEW

(Compatible with BLX Implants)

Regular Base (RB)		30410-01	30410-02	30410-03	30410-04	30410-05	30410-06
Wide Base (WB)		30410-01	30410-02	30410-03	30410-04	30410-05	30410-06

TISSUE LEVEL (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

3.5mm Narrow Neck Crossfit (NNC)		30405-01	30405-02	30405-03	30405-04	30405-05	30405-06
4.8mm Regular Neck (RN)		30406-01	30406-02	30406-03	30406-04	30406-05	30406-06
6.5mm Wide Neck (WN)		30407-01 (1.5mm)	30407-02	30407-03	30407-04	30407-05	30407-06

ZIMMER

INTERNAL HEX CONNECTION (Compatible with Tapered Screw-Vent®, Screw-Vent®, Trabecular Metal™ Implants and AdVent® 5.7mm Platform Implant)

3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06

SPLINE CONNECTION (Compatible with Spline Reliance and Cylinder Implants)

3.25mm Diameter Platform		30205-01	30205-02	30205-03	30205-04	30205-05	30205-06
3.75mm Diameter Platform		30206-01	30206-02	30206-03	30206-04	30206-05	30206-06
4.0mm Diameter Platform		30206-01	30206-02	30206-03	30206-04	30206-05	30206-06
5.0mm Diameter Platform		30207-01	30207-02	30207-03	30207-04	30207-05	30207-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Accessory Components

PART #	DESCRIPTION		PRICE
PROCESSING COMPONENTS			
30001-01	Zero Retention Insert (Gray), (QTY. 4)		\$25.00
30002-01	Low Retention Insert (Blue), (QTY. 4)		\$25.00
30003-01	Medium Retention Insert (Pink), (QTY. 4)		\$25.00
30004-01	High Retention Insert (Clear), (QTY. 4)		\$25.00
30027-01	Low Retention Limited Range Insert (Aqua), (QTY. 4)		\$25.00
30028-01	Medium Retention Limited Range Insert (Magenta), (QTY. 4)		\$25.00
30029-01	High Retention Limited Range Insert (Gold), (QTY. 4)		\$25.00
30012-01	Processing Insert (QTY. 4)		\$25.00
08514	Block-Out Spacer (QTY. 20)		\$20.00
30013-01	Denture Attachment Processing Assembly (QTY. 4)		\$60.00
30014-01	3.35mm Abutment Analog (QTY. 4)		\$75.00
30015-01	4mm Abutment Analog (QTY. 4)		\$75.00
30016-01	5mm Abutment Analog (QTY. 4)		\$75.00
30017-01	Impression Coping (QTY. 4)		\$75.00
30018-01	Processing Spacer (QTY. 4)		\$34.00

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM


Accessory Components (Continued)

TOOLS		
30021-01	LOCATOR R-Tx Retention Insert Tool (Dual Sided) 	\$55.00
08007	.050"/1.25mm Hex Driver with thumb knob (Shaft 6.4mm) 	\$86.00
08008	.050"/1.25mm Hex, Latchlock Driver (Shaft 12.7mm) 	\$53.00
04936	.050"/1.25mm Hex, Torque Wrench Insert (Shaft 6.4mm) 	\$53.00
PART #	DESCRIPTION	PRICE
PATIENT EDUCATION		
09665	LOCATOR R-Tx Patient Education Model	\$180.00
30053-01	LOCATOR R-Tx 4x Demonstration Model	\$70.00

09665



LOCATOR R-Tx®
REMOVABLE ATTACHMENT SYSTEM



VISIT **ZESTDENT.COM**
24 HOURS A DAY | 7 DAYS A WEEK

TO PLACE ORDERS ONLINE OR TO FIND
THE LATEST SUPPORT TOOLS FOR
LOCATOR R-Tx®

**ZD ZEST DENTAL
SOLUTIONS®**

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM

DISTRIBUTED BY: