



Have You *Heard* The News?


Now Available For STRAUMANN® Implants!




LOCATOR F-Tx
FIXED ATTACHMENT SYSTEM




LOCATOR®




LOCATOR R-Tx
REMOVABLE ATTACHMENT SYSTEM

Genuine **LOCATOR**® Attachment Systems Now Available For:

STRAUMANN® IMPLANTS

TISSUE LEVEL • NNC • RN • WN
BONE LEVEL • NC • RC


From The Original Manufacturer **ZEST DENTAL SOLUTIONS**®
And The Zest Distribution Network!

U.S. Retail Product Catalog and Price List for **Straumann Connections**

**ZD ZEST DENTAL
SOLUTIONS**®

ZEST | DANVILLE MATERIALS | PERIOSCOPY

LOCATOR® LOCATOR R-Tx® AND LOCATOR F-Tx® ATTACHMENT SYSTEMS FOR STRAUMANN® IMPLANTS *Abutment Ordering Information*



TO ORDER, YOU WILL NEED THE FOLLOWING INFORMATION:

- 1.** Choose Straumann Implant Type
- 2.** Implant diameter
- 3.** Connection type and platform diameter
- 4.** Cuff height selection criteria

LOCATOR ATTACHMENT SYSTEM: \$119.50 (Abutment Only. Male Processing Pack of 2 Sold Separately \$71.00)

DESCRIPTION	CUFF HEIGHT					
	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
TISSUE LEVEL CONNECTION						
NNC	02056	02057	02058	02059	02060	02061
RN	08742	08743	08744	08745	08751	08752
WN	08870	08871	08872	08873	08874	
BONE LEVEL CONNECTION						
NC		08422	08423	08424	08425	08426
RC	08170	08171	08172	08173	08174	08175

LOCATOR R-Tx ATTACHMENT SYSTEM: \$155 (Includes Abutment and Processing Components)

DESCRIPTION	CUFF HEIGHT					
	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
TISSUE LEVEL CONNECTION						
NNC	30405-01	30405-02	30405-03	30405-04	30405-05	30405-06
RN	30406-01	30406-02	30406-03	30406-04	30406-05	30406-06
WN	30407-01 (1.5MM)	30407-02	30407-03	30407-04	30407-05	30407-06
BONE LEVEL CONNECTION						
NC	30401-01	30401-02	30401-03	30401-04	30401-05	30401-06
RC	30402-01	30402-02	30402-03	30402-04	30402-05	30402-06

LOCATOR F-Tx ATTACHMENT SYSTEM: \$249 (Includes Abutment and Processing Components)

DESCRIPTION	CUFF HEIGHT					
	1.5mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
TISSUE LEVEL CONNECTION						
NNC		10405-02	10405-03	10405-04	10405-05	10405-06
RN		10406-02	10406-03	10406-04	10406-05	10406-06
WN		10407-02 (2.5MM)	10407-03	10407-04	10407-05	10407-06
BONE LEVEL CONNECTION						
NC	10401-01	10401-02	10401-03	10401-04	10401-05	10401-06
RC	10402-01	10402-02	10402-03	10402-04	10402-05	10402-06